

# Overcoming Triple Quadrupole Selectivity Challenges by Use of High-Resolution Accurate Mass Spectrometry to Enable Sensitive Quantitation of a Small Molecule Therapeutic in Rabbit Plasma

Daniel Mulvana<sup>1</sup>, Robert Sturm<sup>1</sup>, Chet Bowen<sup>2</sup>, John Buckholz<sup>1</sup>, Raymond Biondolillo<sup>1</sup>, Kathlyn McCardle<sup>1</sup>, Chris Evans<sup>2</sup>, Barry Jones<sup>1</sup>

<sup>1</sup>Q<sup>2</sup> Solutions, 19 Brown Road, Ithaca, NY, USA. <sup>2</sup>GlaxoSmithKline, King of Prussia, PA, USA.



a Quintiles Quest Joint Venture

## Overview

### Purpose:

- Validated triple quadrupole rat and human plasma methods were used as starting points for the development of a rabbit plasma method for the small molecule GSKA. The lower limit of quantitation (LLOQ) for the rat plasma method was 250 pg/mL using 25 µL of plasma and the LLOQ for the human plasma method was 10 pg/mL using 200 µL of plasma.
- A lower LLOQ (1 pg/mL) was desired for the rabbit plasma method; however, interfering chemical noise due to isobaric interference hindered this effort.
- High-resolution accurate mass spectrometry (HRAMS) was investigated as an alternative technique to enable higher sensitivity by increased selectivity.

### Method:

- GSKA was extracted from rabbit plasma using a double liquid-liquid extraction.
- Chromatographic separation was performed by ultra-high-performance liquid chromatography (UHPLC).
- Heated electrospray ionization (HESI) was used for detection on a Q Exactive (Thermo Scientific) mass spectrometer. Higher energy collision-induced dissociation (HCD) was used to enable tandem mass spectrometric (MS/MS) detection.


### Results:

- An LLOQ of 2 pg/mL using 200 µL of rabbit plasma was achieved.
- The method was validated to full regulatory compliance.

## Introduction

- Triple quadrupole rat and human plasma methods for the small molecule GSKA were challenged by high chemical noise due to isobaric interference. Efforts to generate cleaner extracts included double liquid-liquid extraction with methyl *tert*-butyl ether (MTBE) and pre-rinsing of the 96-well plate with extraction solvent.


### Triple Quadrupole Selected Reaction Monitoring (SRM) Transition Noise (1 pg/mL)



- HRAMS using a Q Exactive was investigated as an alternative technique. Initial efforts using MS1 multiplexed targeted selected ion monitoring (MSX) were unsuccessful due to the isobaric interference, even at 140,000 resolving power.
- Improved selectivity, and thereby sensitivity, was achieved using HCD to enable MS/MS detection.
- Using HCD and 80 µL of rabbit plasma, a 10 pg/mL LLOQ was achieved. Additional method refinement enabled a LLOQ of 2 pg/mL using 200 µL of rabbit plasma.

	Exact Mass	[M+H] <sup>+</sup>
GSKA	390.065 Da	<i>m/z</i> 391.073
SIL IS	397.088 Da	<i>m/z</i> 398.096

### Targeted Multiplexed MSX (MS1) vs HCD (MS2) on the Q Exactive (2.5 ng/mL)



- High baseline noise from an isobaric interferent was observed in the analyte extracted-ion chromatogram (XIC) channel even when a resolving power of 140,000 was used.
- Negligible baseline noise was observed in the analyte channel when HRAMS selectivity was applied to HCD generated product ion spectra at 17,500 resolving power.
- Note: RT differed from the final method because an Acquity UPLC HSS T3 1.8 µm, 2.1 x 50 mm column with an analytical flow rate of 500 µL/min was used.

## Method Overview

### Sample Extraction

- The triple quadrupole method was originally validated in rat plasma using a hydrophilic-lipophilic-balanced (HLB) solid-phase extraction (LLOQ = 0.25 ng/mL). Background interferences prevented reduction of the LLOQ to levels needed for rabbit plasma. The method was modified to use double liquid-liquid extraction with MTBE to improve recovery.

### Scan Rate vs Chromatography

- Orbitrap transient times as a function of required resolving power are displayed below.
- The validated UHPLC triple quadrupole method exhibited a baseline peak width of ~3 sec.
- At low analyte concentrations on the Q Exactive, the baseline peak width was decreased due to limits of detection at the shoulders of the peak elution profile.
- To facilitate adequate peak sampling on the Q Exactive, chromatography was modified to generate a baseline peak width of ~6 sec for mid/high analyte concentrations by reducing flow rate and changing column particle size from 1.7 to 2.5 µm.

### Q Exactive Transient Times for Various Resolution Settings

Resolving Power at <i>m/z</i> 200	Approximate Scan Speed (Hz)	Approximate Scan Time (ms)	Transient Length (ms)	Suggested Maximum Injection Time (ms)
17,500	13	77	64	50
35,000	7	145	128	110
70,000	3	290	256	240
140,000	1.5	580	512	500

- The approximate scan time was matched with the best resolution for analysis based on obtaining at least 10 scans across the chromatographic peak.
- The maximum injection time was matched to the orbitrap transient length to maximize sensitivity and minimize duty cycle.

## Method Development Challenges

### Mobile Phase Considerations

- The validated rat plasma method used a basic mobile phase (0.1% ammonium hydroxide/acetonitrile). Improved sensitivity was achieved using 0.5% ammonium hydroxide/methanol.
- Sodium adducts were observed to contribute to reduced analyte response.
- Previous in-house data showed methanol solvent in glass reservoirs contains high levels of sodium. We now routinely use mobile phase solvents shipped in NOWPak (HDPE/PTFE lined).
- Additionally, use of Teflon mobile phase reservoirs reduces sodium adduction. These are used for some sensitive methods in-house.

### Mass Spectrometric Detection

- The Q Exactive enabled very high resolution mass detection, up to 140,000 resolving power. This is usually sufficient to resolve most isobaric interferences seen under unit mass resolution conditions.
- Targeted MSX was evaluated but was not found to sufficiently reduce background interference using 3 millimass unit (mmu) processing.
- HCD (MS2) fragmentation provided added selectivity, enabling low detection limits.

## Method

### Sample Preparation

- Wash a 2-mL round 96-well plate by adding 500 µL of MTBE to remove possible contaminants. Dry plate.
- Add plasma sample (200 µL) plus SIL IS to the 96-well plate.
- Add 800 µL of MTBE to each well using the Microlab STAR (Hamilton Company).
- Mix on hematology mixer for 10 min.
- Transfer 600 µL of the organic layer to a clean collection plate and evaporate.
- Repeat Steps 3 and 4 and transfer to same evaporation plate.
- After evaporation, reconstitute in water/methanol.

### Chromatographic Conditions

LC System: Dionex/UltiMate 3000 RSLC (Thermo Scientific)  
 Column: XBridge BEH C<sub>18</sub> (2.1 x 50 mm; 2.5 µm particle size; Waters)  
 Column Temperature: 55 °C  
 Injection Volume: 30 µL  
 Mobile Phases\*:  
 A: 1000:5 Water (MilliQ)/Ammonium Hydroxide  
 B: Methanol (NOWPak)

Flow Rate: 250 µL/min  
 Retention Time: 2.80 min  
 \*Teflon mobile phase bottles were used in place of glass mobile phase bottles.  
 Gradient:

Time (min)	Flow Rate (µL/min)	Mobile Phase B (%)
0.00	250	45
3.00	250	95
4.00	250	95
5.00	250	45
6.50	250	45


### Mass Spectrometry Conditions

Instrument: Q Exactive  
 Source: HESI-II  
 Polarity: Positive  
 Spray Voltage: 3.8 kV  
 Sweep Cone Pressure: 0 psi  
 Auxiliary Gas Pressure: 20 psi  
 S-Lens Radio Frequency Level: 50%  
 Capillary Temperature: 370 °C  
 Heater Temperature: 400 °C

### Acquisition Mode

	Targeted MSX (MS1)	HCD (MS2)
Quadrupole Isolation Width	±0.5 <i>m/z</i> units	±2.0 <i>m/z</i> units
Automatic Gain Control Target	5 x 10 <sup>5</sup>	1 x 10 <sup>6</sup>
Maximum Injection Time	240 ms	50 ms
Normalized Collision Energy	-	40%
Orbitrap Resolution	140,000	17,500
Accurate Mass Processing (±3 mmu mass tolerance)	Analyte: <i>m/z</i> 391.073 SIL IS: <i>m/z</i> 398.096	Analyte: <i>m/z</i> 158.976 SIL IS: <i>m/z</i> 158.976


### Is the SRM Interference Observed and Resolved in the Q Exactive HCD Data?



- HRAMS enabled the analyte and SIL IS to be separated from the interferent ion to improve selectivity and overall sensitivity.
- The difference between the HCD product ion and interferent ion was only 0.141 Da.
- Unit resolution mass spectrometers were unable to resolve the interferent from the analyte and SIL IS ions.


## Method Development Changes to Improve Sensitivity

### Mobile Phase: Changing Mobile Phase B to 100% Methanol




- A 2-fold gain in sensitivity was observed by changing Mobile Phase B to 100% methanol.
- Note: RT was different than the final assay because an Acquity UPLC HSS T3 1.8 µm, 2.1 x 50 mm column with an analytical flow rate of 500 µL/min was used.

### Mobile Phase: Increasing Percentage of Ammonium Hydroxide in Mobile Phase A







- A 1.5-fold gain in sensitivity was observed by increasing the ammonium hydroxide concentration in Mobile Phase A.
- Note: RT was different than the final assay because a Waters XBridge BEH C<sub>18</sub> 2.5 µm, 2.1 x 50 mm column with an analytical flow rate of 100 µL/min was used.

### Glass Bottles vs Teflon Bottles as Mobile Phase Reservoirs




- Sodium adduction was decreased by use of Teflon mobile phase reservoirs for mobile phase solutions.
- A 2-fold gain in sensitivity was observed in UHPLC runs on the Q Exactive.
- Note: These data were collected on a TSQ Vantage MS equipped with the same HESI-II source in Q1 only mode.

## Results



Excellent selectivity and negligible carryover were observed with this method on the Q Exactive with HCD.



Excellent linear dynamic range was observed with this method on the Q Exactive with HCD.

### Inter-Day Method Accuracy and Precision

	GSKA Concentration (pg/mL)				
	LLOQ QC 2.00	QC1 6.00	QC2 60.0	QC3 800	QC4 1600
Mean	2.29	6.41	62.8	799	1570
SD	0.178	0.325	1.88	26.3	61.5
CV (%)	7.8	5.1	3.0	3.3	3.9
RE (%)	14.5	6.8	4.7	-0.1	-1.9
n	18	18	18	18	18

- Quality control (QC) sample statistics indicate good accuracy and precision from three method validation runs across the linear dynamic range of 2.00 to 1600 pg/mL.
- The LLOQ in rabbit plasma was established at 2 pg/mL because adequate precision was not demonstrated at 1 pg/mL (data not shown).

### Inter-Lot Accuracy and Precision

	GSKA Concentration (pg/mL)						Inter-lot
	Plasma Lot	Lot 1	Lot 2	Lot 3	Lot 5	Lot 6*	
LLOQ 2.00 pg/mL	Replicate 1	2.01	1.80	2.11	1.94	2.55	1.65
	Replicate 2	2.06	1.66	1.86	1.46	2.12	1.69
	Replicate 3	1.79	2.33	2.22	2.13	2.00	1.57
	Mean	1.95	1.93	2.06	1.84	2.22	1.64
SD	0.144	0.353	0.184	0.345	0.289	0.0611	
CV (%)	7.4	18.3	8.9	18.8	13.0	3.7	
RE (%)	-2.5	-3.5	3.0	-8.0	11.0	-18.0	
n	3	3	3	3	3	3	
ULOQ 2000 pg/mL	Replicate 1	2290	1960	1930	1900	1930	1980
	Replicate 2	2290	1930	1880	1970	1890	1970
	Replicate 3	2370	1960	1900	1900	1900	1920
	Mean	2320	1950	1900	1950	1910	1960
SD	46.2	17.3	25.2	47.3	20.8	32.1	
CV (%)	2.0	0.9	1.3	2.4	1.1	1.6	
RE (%)	16.0	-2.5	-5.0	-2.5	-4.5	-2.0	
n	3	3	3	3	3	3	

- \*Plasma lot with 2% hemolysis.
- Adequate inter-lot accuracy and precision and inter-lot selectivity were observed when processing the HCD data with a ±0.003 Da XIC mass tolerance.
- QC sample statistics indicate good accuracy and precision from six individual lots of rabbit plasma, including one lot with 2% hemolysis.
- In addition, there were no chromatographic interferences observed among the same six individual lots of rabbit plasma, including one lot with 2% hemolysis.

## Conclusions

- Sensitivity was a challenge in the development of the GSKA method in rabbit plasma due to isobaric interference and sodium adduction.
- Isobaric interference was not resolvable by MS1 high-resolution mass spectrometry alone. HCD (MS2) fragmentation was required for selectivity and sensitivity.
- Sodium adduction was minimized by use of non-glass mobile phase solvent reservoirs.
- The method for GSKA in rabbit plasma was successfully developed, validated, and used in regulated preclinical sample analysis.

Acknowledgment: Suzanne Spencer for poster creation